# Instructions for Use

## Sterile Metal Bone Screw

Metal Bone Screws are used either to fixation plates on bones or as big screws to hold bone fragments together. The screws are differentiated by the manner in which they are inserted into bones, their function, their size, and the type of bone they are intended for.

**I. Material**  
Metal Bone Screws are made of Titanium alloy.

**II. Intended Purpose**  
Metal Bone Screw is intended to be used alone or in combination with bone plates for temporary fixation, correction and stabilization of limb fractures.

**III. Indication**  
Metal Bone Screw is indicated for internal fixation of limb fractures.

**IV. Contraindications**  
(1) General or local infection, osteomyelitis.  
(2) Patient can't tolerate operation or anesthesia because of poor health.  
(3) Mental disease, systemic neurological disease,  
(4) Pregnancy, metabolic disorders caused by calcified tissue.  
(5) Division of the tissue can lead to infection or skin necrosis due to poor local soft tissue conditions.  
(6) Severe osteoporosis prevents effective fixation.  
(7) Patients with influence factors for fracture healing (e.g., diabetes, infection, etc.).  
(8) Patients with high degree of comminution, III degree open, severe bone defect and severe damage to surrounding soft tissue.

**V. Intended user**  
The operation should be performed by professional trained surgeons according to the IFU and surgical technique guide strictly. The surgeon should have a thorough understanding of the intended purpose of the product and the surgical techniques used, which should be trained and certified by relevant institutions (e.g., Association for Internal Fixation Research, AO/ASIF).

**VI. Patient Group**  
Adults who have the surgical indications and who have no surgical contraindications and can tolerate surgery.

**VII. Clinical benefits**  
Metal Bone Screw is used for internal fixation of bone fracture to achieve bone union with expected union rate no less than 95%.

**VIII. Potential Side Effects**  
1. Implant displacement screw penetration, loosening or implant fracture.  
2. Nonunion, delayed union or malunion or bony reaction.  
3. Adverse tissue reactions, allergy/hypersensitivity reaction.  
4. Infection.  
5. Soft tissue damage or neuro-vascular damage.  
6. Hematoma reactions, osteolysis or bone necrosis.  
7. Pain.  
8. Impingement or poor joint mechanics.  
9. Wear or corrosion of metal implants.

**IX. Warnings and Precautions**  
1. Preoperative planning and adequate training of patients, adequate training, experience in the choice and placement of implants and the decision to leave or remove implants postoperatively, rests with the surgeon.  
2. The surgeon should discuss the expectations of surgery inherent in the use of the product with the patient preoperatively. Particular attention should be paid to a discussion postoperatively and the necessity for periodic re-evaluation of the patient's condition.  
3. Suitable implant should be chosen for specific fracture and implant indication. Choose a suitable metal locking bone screw according to preoperative X-ray film. Select a supporting instrument set and metal locking bone screw according to the chosen locking bone plate. Surgery should be finished with the guidance of Manual of Instructions Recommended by the AO-ASIF Group and other relevant information.  
4. The correct selection of the product is extremely important. The product should be used in the correct anatomic location, consistent with accepted standards for internal fixation. Failure to use the appropriate product for the application may result in a premature clinical failure. Failure to use the proper component to maintain adequate fixation and provide rigid fixation may result in loosening, bending, or fracturing of the product and/or bone.  
5. Careful preoperative planning is a requirement. The surgeons should choose the proper implants according to preoperative images and by using relating measuring instruments such as depth gauges and trials. Instruments which intended for use with the implant should be used in order to obtain a safe combination.  
6. In order to avoid or minimize specific risks associated with implantation, surgeons should have been trained, and the operation should be performed strictly according to the surgical technique guide.  
7. Intraoperative decisions should not be revised. The risk of reuse is as follows:  
(1) Metal strength decreases and fracture occurs easily after being implanted.  
(2) As contaminants and bacteria remaining on the surface, re-implantation would result in infection and rejection.  
(3) It can damage the surface anodic oxide film and affect the biological properties of the product.  
(4) Careful handling and storage of the product is required. Scratching or damage to the component can significantly reduce the strength and fatigue resistance of the product.  
8. Generally, there is no limitation on substance used in the clinical setting, to which the implant might be exposed.  
9. Once applied, the product should never be reused. Although it may appear undamaged, previous stresses may have created imperfections that could reduce its service life.  
10. In extreme angles together with small bending radii must be avoided because of the potential risk for fracture.  
11. Excessively aggressive use of bending instruments can cause regardless macroscopic damage to the implant (indentations, engraved screw holes, etc.). In such cases, the implant must be exchanged for a new one.  
12. Unless specifically indicated in the surgical technique guide, the modification of size, shape or surface texture is prohibited.

**14. Metal Bone Screw has not been evaluated for safety in the MR environment. It has not been tested for heating or unwarmed movement in the MR environment. The safety of Metal Bone Screw in the MR environment is unknown. Performing an MR exam on an implant who has this medical device may result in injury or device failure.**

15. Surgeons should verify whether the patients are positioned correctly using adequate imaging technology.  
16. Patients should be advised to seek medical opinion before entering potentially adverse environments that could affect the performance of the implant.  
17. Patients should avoid physical activity or load bearing after surgery is prohibited.  
18. The patient should be advised to report any unusual changes of the operated site to his surgeon. The patient should be closely monitored if a change at the fixation site has been detected. The surgeon should evaluate the possibility of subsequent clinical failure, and discuss with the patient the need for any measures deemed necessary to avoid re-healing.  
19. Postponed or interrupted use may bring about implant loosening or crack; therefore, after bone healing, it's suggested to take out the implant if there is implant crack risk, which may bring medical accident.  
20. Considerations for removal of the implant after healing. If the device is not removed after the completion of its intended use, any undesired complications may occur (1) corrosion, with localized tissue reaction or pain; (2) migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) pain, discomfort, or abnormal sensations due to the presence of the device; (6) possible increased risk of infection; and (7) bone healing delay. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management and wound care. If the patient is older and has a low activity level, the surgeon may choose not to remove the implant unless the risks involved with a second surgery.

**X. Sterility**  
The implant is provided as sterile. It is sterilized with  $\gamma$  ray, and the validity period is 5 years.

	Manufacturer
	Use-by date
	Do not re-use
	Do not use if package is damaged
	Caution
	Keep dry
	Importer
	Double sterile barrier system

**XII. Storage**  
The sterilized product should be stored in a clean environment, protected from direct sunlight and pests.

**XIII. Precautions**  
Avoid collision and compression during storage and transportation.

**XIV. Use of Original Products**  
Implants and instruments are designed to be used together. The use of products from other manufacturers along with Double Medical products may involve product risks, and corrosion of the material and misalignment of implant and instruments, impeding functionality, thereby endangering the patient, user or third parties.

**XV. Handling Information**  
Metal Bone Screws are made of Titanium alloy. The material is biocompatible as widely used in the industry, corrosion-resistant and non-toxic in the biological environment. It produces negligible artifacts by X-ray and CT.

**XVI. Surgical Technique**  
(1) Take a X-ray picture of the fractured bone, including adjacent joints.  
(2) Determine the type of the original bone status and the correction of bone malformation are the main purpose of the treatment.  
(3) Assess the condition of soft tissue and exam the neurological and vascular functions.  
(4) Surgical Approach  
The length of incision is decided as needed. Cut the skin, subcutaneous tissue, fascia and peritosteum. Expose the fracture and remove to protect nerve and blood vessel.  
(5) Determination  
Determine plate length and bend plate if plate is used.  
(6) Drilling  
Conduct a pre-drill with an invasive predrill or orthopaedic traction table. Choose patient position according to fracture type.  
(7) Restoration and Temporary Fixation  
Use a reduction pin and Kirschner wire to restore. Restoration result can be examined when needed.  
(8) Placement of Implant  
Place a nail on the surface of fractured bone and select a bone plate with proper length and shape.  
(9) Drilling  
Use a hand drill using a proper drill and guide and compress as needed.  
(10) Inserting Screw Length  
Use a detector to detect the depth of holes and the depth is the length of screws. When detecting, the detector should reach the cortex on the opposite side and hook it.  
(11) Tapping  
Use a tapping rod with screw tap tip to tapping.  
(12) Screw Insertion  
Before inserting the first screw, ensure the fracture has been restored. Tap the screw in. If the screw is long or the cortex is too thick, normal saline can be applied to lower the temperature.  
(13) Tightening  
Expose the plate and screws from the original incision and strip tissues on the surface of the bone plate and screws. Use a hexagonal socket screw driver to unscrew, remove all the screws and finally remove the plate.  
(14) Final Use  
Finally connect the hexagonal socket screw driver and the screw to make them in line in case of gear slip.

**The handling of gear slip or screw crack**  
(1) When gear slip occurs and it is impossible for the hexagonal socket screw driver to remove the screw, connect a slip extractor to an electrical bone drill and insert vertically into the hexagonal socket and the screw. Tighten the screw and rotate the gear slip extractor.  
(2) When screw crack occurs, connect a crack extractor to an electrical bone drill and drill the bone cortex around the screw antiockwise and the screw can be removed.

**14. Disposal of Implant**  
Implants removed from any body should be handled in an anti-pollution way according to hospital protocol so as to prevent cross infection.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

**XVII. Manufacturer Contact Information**  
Manufacturer: DouMe Medical Technology Inc.  
Legal Address: No. 18, Shantianbin East Road, Haicang District, 361026, Xiamen, Fujian, PEOPLES REPUBLIC OF CHINA  
Telephone Number: +86 592 6087101  
SSCP (Summary of safety and clinical performance): Eudamed (https://ec.europa.eu/tools/eudamed/screen/home)  
e-IFU: <https://www.doublemedical.com>

Note: If any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Notified Body: TÜV SÜD Product Service GmbH  
Address: Riederstrasse 65, 80339 Munich, Germany

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Fabricante: DouMe Medical Technology Inc.
Organismo notificado: TÜV SÜD Product Service GmbH
Shanghai International Holding Corp. GmbH (Europa)