



The latest solution for the fusion of tibiotalar joint.



Combi-holes design is convenient for compression and locking.



Two screws on the shaft can be inserted from tibia to talus which enhances the fixation.

Indications

– Fusion of tibiotalar joint. The plate is recommended to place on the anterior lateral tibia.



Warning


This flyer is just for understanding the specific product features. For clinical usage, please refer to the surgical guide. Instruction by experienced surgeon is highly recommended.



3.5mm Locking Screw, Self-drilling



3.5mm Locking Screw, Self-tapping



3.5mm Cortex Screw



4.0mm Cancellous Screw, Half Threaded